

# R documentation

of 'F:/documents/these/stat\_R/R' etc.

August 5, 2015

## R topics documented:

BrainMetabolism-package . . . . .	1
auc . . . . .	2
calibration . . . . .	2
CMRGluc.calc . . . . .	3
CMRO2.calc . . . . .	4
correction.Temp . . . . .	5
correction.TPO2 . . . . .	5
IO2.calc . . . . .	6
L.calc . . . . .	7
noise.na . . . . .	8
OGI.calc . . . . .	8
polyfit . . . . .	9
polyval . . . . .	9
roll.funct . . . . .	10

---

BrainMetabolism-package

*calculate brain metabolism rate from extracellular concentrations*

---

## Description

calculate CMRO2, CMRgluc, mitochondrialPO2 ...  
provide functions to process data from biosensors

## Details

Package: BrainMetabolism  
Type: Package  
Version: 1.0  
Date: 2015-03-03  
License: MIT

**Author(s)**

Baptiste Balanca  
 Lyon Neuroscience Research Center  
 Team TIGER  
 Maintainer: baptiste balanca <baptiste.balanca@gmail.com>

**References**

Piilgaard Lauritzen JCBFM (2009) 29, 1517  
 Gjedde et al JCBFM (2005) 25(9), 1183  
 Gjedde et al JCBFM (2000) 20(4), 747  
 Du et al JCBFM 2012 32(9)  
 Gruetter et al J Neurochem 1998 70(1)  
 balanca et al (2016) in preparation

**Examples**

```
MyTPO2<-23 #mmHg
MyLDF<-95 #percent
MyL<-L.calc(CMR=219, TPO2=MyTPO2, P50=36, h=2.7, Ca=8, cbf=53 )
MyCMRO2<-CMRO2.calc(LDF=MyLDF, MyTPO2, P50=36, h=2.7, Ca=8, L=MyL, cbfbase=53)

MyGlucBrain<-1 #mM
MyGlucPlasma<-6 #mM
MyCRMgluc<-CRMgluc.calc(MyGlucbrain, Vd=0.77, Kt=13.4, Tmax=1.35, Gplasma=MyGlucPlasma)
```

---

auc	<i>area under the curve</i>
-----	-----------------------------

---

**Description**

take a numeric vector and return the area under the curve (AUC)

**Usage**

```
auc(data, pos = TRUE)
```

**Arguments**

data	: a numeric vector
pos	: a logical value. if TRUE (default) AUC is calculated on positive values, if FALSE on negative values.

**Value**

auc : area under the curve  
 max/min : maxiaml value (or minimal if pos=FALSE)

---

calibration	<i>Sensor calibration</i>
-------------	---------------------------

---

**Description**

Fit an  $n^{th}$  degree polynom

$$(y = C_n X^n + C_{n-1} X^{n-1} + \dots + C_1 X + C_0)$$

to vectors x=volt and y=concentration

**Usage**

```
calibration(volt, mol, order = 1)
```

**Arguments**

volt : a numeric vector of the biosensor voltage

mol : a numeric vector of the molecule concentration in the medium (x and y must have the same length)

order : a numeric value for the polynomial degree. default is one.

**Value**

Coef : coefficients in ascending order (i.e. C0, C1, C2, ..., Cn )

R2 : goodness of fit

---

CMRGluc.calc	<i>CMRGluc calculation</i>
--------------	----------------------------

---

**Description**

take extracellular glucose concentration and return brain metabolic rate of glucose, using a reversible Michaelis-Menten model equation:

**Usage**

```
CMRGluc.calc(Gbrain, Vd = 0.77, Kt = 13.4, Tmax = 1.35, Gplasma = 6)
```

**Arguments**

Gbrain : a vector/numeric value of extracellular glucose concentration in mmol/L, time decay

Vd : glucose brain space diffusion (default is 0.77 ml/g)

Kt : glucose apparent maximal transport rate (default is 13.4 mmol/L)

Tmax : the apparent maximal transport rate (default is 1.35 micromol/g/min)

Gplasma : plasma glucose concentration (default=6mmol/L)

## Details

$$G_{brain} = V_d \frac{\left( \frac{T_{max}}{CMR_{gluc}} - 1 \right) \times G_{plasma} - K_t}{\frac{T_{max}}{CMR_{gluc}} + 1}$$

## Value

CMRGlucose : numeric value of Cerebral metabolic rate of glucose in micromol/g/min

## References

Du et al JCBFM 2012 32(9)

Gruetter et al J Neurochem 1998 70(1)

---

CMRO2.calc

*CMRO2 calculation*

---

## Description

take tissue oxygen pressure (tPO2) and cerebral blood flow (relative value, e.g. laser doopler lowmetry BPU)

## Usage

```
CMRO2.calc(LDF, TPO2, P50 = 36, h = 2.7, Ca = 8, L = 4.03,
           cbfbase = 53)
```

## Arguments

LDF	: a vector/numeric value of LDF (percentage of baseline, i.e: basal value is 100 % )
TPO2	: a vector/numeric value of brain oxyngen pressure (mmHg)
P50	: half-saturation tension of hemoglobine (default is 36 mmHg)
h	: hill's coeficient
Ca	: oxygen arterial concentration (default is 8 micromol/ml)
L	: effective diffusion coefficient of oxygen in brain tissue, default is 4.03 micromol/100g/mmHg, but one should use the L.calc function to calculate it from their data.
cbfbase	: basal expected value of CBF (default is 53 ml/100g/min wich was used to calculate L) LBF and TPO2 must be the same length

## Details

$$PbtO_2 = P_{50} \cdot \sqrt[h]{\frac{2 \cdot C_a \cdot CBF}{CMRO_2} - 1} - \frac{CMRO_2}{2 \cdot L}$$

**Value**

CMRO2 : vector/numeric value of Cerebral Metabolic Rate of Oxygen in micromol/100g/min

**References**

Gjedde et al JCBFM (2005) 25(9), 1183

Piilgaard et al JCBFM (2009) 29, 1517

---

correction.Temp	<i>Temperature correction</i>
-----------------	-------------------------------

---

**Description**

Biosensor enzymatic reaction, that underpin amperometric measures, has a sigmoid relation to temperature :

$$m(x, P) = \frac{P_1 + (P_2 - P_1)}{(1 + \exp((P_3 - x)/P_4))}$$

parameters are different for each enzyme and has been measured in vitro.

This function correct biosensor signal depending on temperature conditions during calibration and experimentation

**Usage**

```
correction.Temp(x, enz, temp.calib = 25, temp.exp = 37)
```

**Arguments**

x	: a numeric vector of the biosensor voltage
enz	: enzyme on the biosensor i.e. "glucose", "lactate", "glutamate", "daao"
temp.calib	: temperature of the medium where sensor has been calibrated the default is 25°C (i.e. room temperature).
temp.exp	: temperature of the medium during experiment the default is 37°C (i.e. animal central temperature)

**Value**

volt.temp.cor : a vector of corected x values for temperature

**References**

balanca et al 2015

---

correction.TPO2	<i>Oxygen Tension correction</i>
-----------------	----------------------------------

---

### Description

Biosensor enzymatic reaction, that underpin amperometric measures, has an asymptotic relation to oxygen tension in the medium :

$$m(PO_2, P) = P_1 + (P_2 - P_1) \times \exp(-\exp(P_3)PO_2)$$

parameters are different for each enzyme and has been measured in vitro

This function correct biosensor signal depending on PO2 conditions during calibration and experimentation

### Usage

```
correction.TPO2(x, enz, TPO2 = 28)
```

### Arguments

x	: a numeric vector of the biosensor voltage
enz	: enzyme on the biosensor i.e. "glucose", "lactate", "glutamate", "daao"
TPO2	: oxygene tension in the medium during the experiment. the default is 30mmHg measured in anesthetized rat brain.

### Value

volt.O2.cor : a vector of corected x values for TPO2

### References

balanca et al 2015

---

IO2.calc	<i>Oxidative index calculation</i>
----------	------------------------------------

---

### Description

take CMRO2 and LDF to give an oxidative index

### Usage

```
IO2.calc(CMRO2, LDF, cbfbasal = 53)
```

### Arguments

CMRO2	: vector/numeric value of Cerebral metabolic rate of oxygen in micromol/100g/min
LDF	: a vector/numeric value of LDF (percentage from baseline, baseline is 100 %)
cbfbasal	: basale expected value of CBF from the litterature (default is 53 ml/100g/min wich was used to calculate L)

**Details**

$$IO2 = CMRO2 / (cbf_{basal} * LDF)$$

**Value**

IO2 : oxidative index, reflect the degree of flow metabolism coupling

**References**

Gjedde et al JCBFM (2000) 20(4), 747

---

L.calc	<i>Calculate the effective diffusion coefficient of oxygen in brain tissue, L.</i>
--------	--

---

**Description**

take CMRO2 and cerebral blood flow (CBF) to calculate the effective diffusion coefficient of oxygen in brain tissue (L)

**Usage**

L.calc (CMR = 219, TPO2, P50 = 36, h = 2.7, Ca = 8, cbf = 53)

**Arguments**

CMR	: Cerebral Metabolic Rate of Oxygen, default is 219 micro mol/100
TPO2	: numeric value of brain oxygen pressure (mmHg)
P50	: half-saturation tension of hemoglobine (default is 36 mmHg)
h	: hill's coefficient
Ca	: oxygen arterial concentration (default is 8 micromol/ml)
cbf	: expected value of CBF (default is 53 ml/100g/min wich was used to calculate L

**Details**

$$PbtO_2 = P_{50} \cdot \sqrt[h]{\frac{2 \cdot C_a \cdot CBF}{CMRO_2} - 1} - \frac{CMRO_2}{2 \cdot L}$$

**Value**

L : numeric value of the effective diffusion coefficient of oxygen in brain tissue micomol/100g/mmHg

**References**

Gjedde et al JCBFM (2005) 25(9), 1183  
 Piilgaard et al JCBFM (2009) 29, 1517

---

<code>noise.na</code>	<i>Remove artifacts from biosensor signal, based on data's standar deviation</i>
-----------------------	--

---

### Description

Remove artifacts from biosensor signal, based on data's standar deviation

### Usage

```
noise.na(data, z = 20, width = 30)
```

### Arguments

<code>data</code>	: a numeric vector
<code>z</code>	: a numeric value. Number of SD over which values should be exculded
<code>width</code>	: size of the window used to roll SD over data (see <code>roll.funct</code> )

### Value

a vector with NA remplacing exculded values

---

<code>OGI.calc</code>	<i>Oxygene Glucose index (OGI)</i>
-----------------------	------------------------------------

---

### Description

take CMRO2 and CMRGlucose to give an OGI

### Usage

```
OGI.calc(CMRO2, CMRGluc)
```

### Arguments

<code>CMRO2</code>	: vector (numeric value) of Cerebral metabolic rate of oxygen in micromol/100g/min
<code>CMRGluc</code>	: a vector (numeric value) of CMRGlucose in micromol/100g/min

### Details

$$OGI = CMRO2/CMRGlucose$$

### Value

OGI



polyfit

*polynomial fit***Description**

Fit a  $n^{th}$  degree polynom

$$(y = C_n X^n + C_{n-1} X^{n-1} + \dots + C_1 X + C_0)$$

to vectors x and y

**Usage**

```
polyfit(x, y, order = 1)
```

**Arguments**

x : a numeric vector  
 y : a numeric vector (x and y must have the same length)  
 order : a numeric value for the polynomial degree. default is one.

**Value**

model formula  
 Coef: coefficients in ascending order (i.e. C0, C1, C2, ..., Cn )  
 R2: Rsquare

polyval

*polynomial evaluation***Description**

evaluation a  $n^{th}$  degree polynom

$$(y = C_n X^n + C_{n-1} X^{n-1} + \dots + C_1 X + C_0)$$

at given values of x

**Usage**

```
polyval(x, coef)
```

**Arguments**

x : a numeric vector  
 coef : a n dimation vector corresponding to the polynom coefficient (ascending order, i.e. C0, C1, C2, ... , Cn )

**Value**

y

---

roll.funct	<i>apply a function FUN on a rooling windows of a vector</i>
------------	--

---

**Description**

apply a function FUN on a rooling windows of a vector

**Usage**

```
roll.funct(data, width, FUN, size = T, ...)
```

**Arguments**

data	: a numeric vector
width	: the size of the rolling window
FUN	: the function to apply
size	: a logical value indicating if the returned vector have the same length as original data. default is TRUE.
...	: additional argument to pass to the FUN

**Value**

a vector with FUN result