

Package ‘BrainMetabolism’

August 23, 2015

Type Package

Title calculate brain metabolism rates from extracellular concentrations

Version 1.0

Date 2015-03-03

Author baptiste balanca

Maintainer baptiste balanca <baptiste.balanca@gmail.com>

Depends R (>= 3.1.0), ggplot2, stats

Imports ggplot2

Description calculate CMRO2, CMRgluc, mitochondrialPO2 ... provide functions to process data from biosensors

License MIT

NeedsCompilation no

R topics documented:

BrainMetabolism-package	2
auc	3
calibration	3
CMRGluc.calc	4
CMRO2.calc	5
correction.Temp	6
correction.TPO2	6
IO2.calc	7
L.calc	8
noise.na	9
OGI.calc	9
polyfit	10
polyval	10
roll.funct	11
Index	12

BrainMetabolism-package

calculate brain metabolism rate from extracellular concentrations

Description

calculate CMRO2, CMRgluc, mitochondrialPO2 ...
provide functions to process data from biosensors

Details

Package: BrainMetabolism
Type: Package
Version: 1.0
Date: 2015-03-03
License: MIT

Author(s)

Baptiste Balanca
Lyon Neuroscience Research Center
Team TIGER
Maintainer: baptiste balanca <baptiste.balanca@gmail.com>

References

Piilgaard Lauritzen JCBFM (2009) 29, 1517
Gjedde et al JCBFM (2005) 25(9), 1183
Gjedde et al JCBFM (2000) 20(4), 747
Du et al JCBFM 2012 32(9)
Gruetter et al J Neurochem 1998 70(1)
balanca et al (2016) in preparation

Examples

```
MyTP02<-23 #mmHg
MyLDF<-95 #percent
MyL<-L.calc(CMR=219, TP02=MyTP02, P50=36, h=2.7, Ca=8, cbf=53 )
MyCMRO2<-CMRO2.calc(LDF=MyLDF, MyTP02, P50=36, h=2.7, Ca=8, L=MyL,cbfbase=53)

MyGlucBrain<-1 #mM
MyGlucPlasma<-6 #mM
MyCRMgluc<-CMRGluc.calc(MyGlucbrain, Vd=0.77, Kt=13.4, Tmax=1.35, Gplasma=MyGlucPlasma)
```

auc	<i>area under the curve</i>
-----	-----------------------------

Description

take a numeric vector and return the area under the curve (AUC)

Usage

```
auc(data, pos = TRUE)
```

Arguments

data : a numeric vector
 pos : a logical value. if TRUE (default) AUC is calculated on positive values, if FALSE on negative values.

Value

auc : area under the curve
 max/min : maxiaml value (or minimal if pos=FALSE)

calibration	<i>Sensor calibration</i>
-------------	---------------------------

Description

Fit an n^{th} degree polynom

$$(y = C_n X^n + C_{n-1} X^{n-1} + \dots + C_1 X + C_0)$$

to vectors x=volt and y=concentration

Usage

```
calibration(volt, mol, order = 1)
```

Arguments

volt : a numeric vector of the biosensor voltage
 mol : a numeric vector of the molecule conentration in the medium (x and y must have the same length)
 order : a numeric value for the polynomial degree. default is one.

Value

Coef : coeficients in ascending order (i.e. C0, C1, C2, ..., Cn)
 R2 : goodness of fit

CMRGluc.calc

CMRGluc calculation

Description

take extracellular glucose concentration and return brain metabolic rate of glucose, using a reversible Michaelis-Menten model equation:

Usage

CMRGluc.calc(Gbrain, Vd = 0.77, Kt = 1.4, Tmax = 1.27, Gplasma = 7.3)

Arguments

Gbrain	: a vector/numeric value of extracellular glucose concentration in mmol/L, time decay
Vd	: glucose brain space diffusion (default is 0.77 ml/g)
Kt	: glucose apparent maximal transport rate (default is 1.4 mmol/L)
Tmax	: the apparent maximal transport rate (default is 1.27 micromol/g/min)
Gplasma	: plasma glucose concentration (default=6mmol/L)

Details

$$G_{brain} = V_d \frac{\left(\frac{T_{max}}{CMR_{gluc}} - 1 \right) \times G_{plasma} - K_t}{\frac{T_{max}}{CMR_{gluc}} + 1}$$

Value

CMRGlucose : numeric value of Cerebral metabolic rate of glucose in micromol/g/min

References

Morgenthaler et al. Neurochem Int 2006, 48
 Ori et al. Anesthesiology 1986, 65(2)
 Duarte et al. Front in Neuroenrgetics 2009, 1
 Gruetter et al J Neurochem 1998, 70(1)

CMRO2.calc

CMRO2 calculation

Description

take tissue oxygen pressure (tPO2) and cerebral blood flow (relative value, e.g. laser doppler lowmetry BPU)

Usage

```
CMRO2.calc(LDF, TPO2, P50 = 36, h = 2.7, Ca = 8, L = 4.03,
           cbfbase = 53)
```

Arguments

LDF	: a vector/numeric value of LDF (percentage of baseline, i.e: basal value is 100 %)
TPO2	: a vector/numeric value of brain oxygen pressure (mmHg)
P50	: half-saturation tension of hemoglobine (default is 36 mmHg)
h	: hill's coefficient
Ca	: oxygen arterial concentration (default is 8 micromol/ml)
L	: effective diffusion coefficient of oxygen in brain tissue, default is 4.03 micromol/100g/mmHg, but one should use the L.calc function to calculate it from their data.
cbfbase	: basal expected value of CBF (default is 53 ml/100g/min wich was used to calculate L) LBF and TPO2 must be the same length

Details

$$PbtO_2 = P_{50} \cdot \sqrt[h]{\frac{2 \cdot C_a \cdot CBF}{CMRO_2} - 1} - \frac{CMRO_2}{2 \cdot L}$$

Value

CMRO2 : vector/numeric value of Cerebral Metabolic Rate of Oxygen in micromol/100g/min

References

Gjedde et al JCBFM (2005) 25(9), 1183
Piilgaard et al JCBFM (2009) 29, 1517

correction.Temp	<i>Temperature correction</i>
-----------------	-------------------------------

Description

Biosensor enzymatic reaction, that underpin amperometric measures, has a sigmoid relation to temperature :

$$m(x, P) = \frac{P_1 + (P_2 - P_1)}{(1 + \exp((P_3 - x)/P_4))}$$

parameters are different for each enzyme and has been measured in vitro.

This function correct biosensor signal depending on temperature conditions during calibration and experimentation

Usage

```
correction.Temp(x, enz, temp.calib = 25, temp.exp = 37)
```

Arguments

x	: a numeric vector of the biosensor voltage
enz	: enzyme on the biosensor i.e. "glucose", "lactate", "glutamate", "daao"
temp.calib	: temperature of the medium where sensor has been calibrated the default is 25°C (i.e. room temperature).
temp.exp	: temperature of the medium during experiment the default is 37°C (i.e. animal central temperature)

Value

volt.temp.cor : a vector of corected x values for temperature

References

balanca et al 2015

correction.TPO2	<i>Oxygen Tension correction</i>
-----------------	----------------------------------

Description

Biosensor enzymatic reaction, that underpin amperometric measures, has an asymptotic relation to oxygen tention in the medium :

$$m(PO_2, P) = P_1 + (P_2 - P_1) \times \exp(-\exp(P_3)PO_2)$$

parameters are different for each enzyme and has been measured in vitro

This function correct biosensor signal depending on PO2 conditions during calibration and experimentation

Usage

```
correction.TPO2(x, enz, TPO2 = 28)
```

Arguments

x : a numeric vector of the biosensor voltage
 enz : enzyme on the biosensor i.e. "glucose", "lactate", "glutamate", "daao"
 TPO2 : oxygene tension in the medium during the experiment. the default is 30mmHg measured in anesthetized rat brain.

Value

volt.O2.cor : a vector of corected x values for TPO2

References

balanca et al 2015

 IO2.calc

Oxidative index calculation

Description

take CMRO2 and LDF to give an oxidative index

Usage

```
IO2.calc(CMRO2, LDF, cbfbasal = 53)
```

Arguments

CMRO2 : vector/numeric value of Cerebral metabolic rate of oxygen in micromol/100g/min
 LDF : a vector/numeric value of LDF (percentage from baseline, baseline is 100 %)
 cbfbasal : basale expected value of CBF from the litterature (default is 53 ml/100g/min wich was used to calculate L)

Details

$$IO2 = CMRO2 / (cbfbasal * LDF)$$

Value

IO2 : oxidative index, reflect the degree of flow metabolism coupling

References

Gjedde et al JCBFM (2000) 20(4), 747

L.calc	<i>Calculate the effective diffusion coefficient of oxygen in brain tissue, L.</i>
--------	--

Description

take CMRO2 and cerebral blood flow (CBF) to calculate the effective diffusion coefficient of oxygen in brain tissue (L)

Usage

L.calc(CMR = 219, TP02, P50 = 36, h = 2.7, Ca = 8, cbf = 53)

Arguments

CMR	: Cerebral Metabolic Rate of Oxygen, default is 219 micro mol/100
TP02	: numeric value of brain oxygen pressure (mmHg)
P50	: half-saturation tension of hemoglobine (default is 36 mmHg)
h	: hill's coefficient
Ca	: oxygen arterial concentration (default is 8 micromol/ml)
cbf	: expected value of CBF (default is 53 ml/100g/min wich was used to calculate L

Details

$$P_{btO_2} = P_{50} \cdot \sqrt[h]{\frac{2 \cdot C_a \cdot CBF}{CMRO_2} - 1} - \frac{CMRO_2}{2 \cdot L}$$

Value

L : numeric value of the effective diffusion coefficient of oxygen in brain tissue micomol/100g/mmHg

References

Gjedde et al JCBFM (2005) 25(9), 1183
 Piilgaard et al JCBFM (2009) 29, 1517

noise.na	<i>Remove artifacts from biosensor signal, based on data's standar deviation</i>
----------	--

Description

Remove artifacts from biosensor signal, based on data's standar deviation

Usage

```
noise.na(data, z = 20, width = 30)
```

Arguments

data	: a numeric vector
z	: a numeric value. Number of SD over which values should be exculded
width	: size of the window used to roll SD over data (see roll.funct)

Value

a vector with NA replacing exculded values

OGI.calc	<i>Oxygene Glucose index (OGI)</i>
----------	------------------------------------

Description

take CMRO2 and CMRGlucose to give an OGI

Usage

```
OGI.calc(CMRO2, CMRGluc)
```

Arguments

CMRO2	: vector (numeric value) of Cerebral metabolic rate of oxygen in micromol/100g/min
CMRGluc	: a vector (numeric value) of CMRGlucose in micromol/100g/min

Details

$$OGI = CMRO2/CMRGlucose$$

Value

OGI

polyfit	<i>polynomial fit</i>
---------	-----------------------

Description

Fit a n^{th} degree polynom

$$(y = C_n X^n + C_{n-1} X^{n-1} + \dots + C_1 X + C_0)$$

to vectors x and y

Usage

polyfit(x, y, order = 1)

Arguments

x : a numeric vector
 y : a numeric vector (x and y must have the same length)
 order : a numeric value for the polynomial degree. default is one.

Value

model formula

Coef: coefficients in ascending order (i.e. C0, C1, C2, ..., Cn)

R2: Rsquare

polyval	<i>polynomial evaluation</i>
---------	------------------------------

Description

evaluation a n^{th} degree polynom

$$(y = C_n X^n + C_{n-1} X^{n-1} + \dots + C_1 X + C_0)$$

at given values of x

Usage

polyval(x, coef)

Arguments

x : a numeric vector
 coef : a n dimation vector corresponding to the polynom coefficient (ascending order, i.e. C0, C1, C2, ... , Cn)

Value

y

roll.funct	<i>apply a function FUN on a rooling windows of a vector</i>
------------	--

Description

apply a function FUN on a rooling windows of a vector

Usage

roll.funct(data, width, FUN, size = T, ...)

Arguments

data	: a numeric vector
width	: the size of the rolling window
FUN	: the function to apply
size	: a logical value indicating if the returned vector have the same length as original data. default is TRUE.
...	: additional argument to pass to the FUN

Value

a vector with FUN result

Index

*Topic **package**

BrainMetabolism-package, [2](#)

auc, [3](#)

BrainMetabolism

(BrainMetabolism-package), [2](#)

BrainMetabolism-package, [2](#)

calibration, [3](#)

CMRGluc.calc, [4](#)

CMR02.calc, [5](#)

correction.Temp, [6](#)

correction.TPO2, [6](#)

I02.calc, [7](#)

L.calc, [8](#)

noise.na, [9](#)

OGI.calc, [9](#)

polyfit, [10](#)

polyval, [10](#)

roll.funct, [11](#)